

IN THE UNITED STATES PATENT AND TRADEMARK OFFICE

Application of:	Rice et al.	Group No.:	1648
Serial No.:	09/576,989	Atty. Docket No.:	56029-4356
Filed:	05/23/2000		
For:	HCV Variants	Examiner:	Wortman, Donna C.

DECLARATION OF DR. KERIL J. BLIGHT UNDER 37 C.F.R. §1.131

I, Dr. Keril J. Blight, declare and state as follows:

1. All of the statements made herein of my own knowledge are true and all statements made on information and belief are believed to be true.
2. I am a co-inventor of U.S. Patent Application No. 09/576,989 for HCV Variants, filed May 23, 2000 (Patent Application).
3. I, along with co-inventor Dr. Charles M. Rice (hereinafter, "we"), conceived of, and reduced to practice the inventions claimed in the Patent Application before March 31, 2000.

4. Specifically, before March 31, 2000 we identified the adaptive mutations that are described in the Patent Application. Those adaptive mutations are referenced in the attached laboratory notebook pages and computer printouts attached as Exhibit A. The terminology used to describe the cell colonies harboring HCV comprising those mutations in Exhibit A (see A25) corresponds to the terminology used in the Patent Application (see Figure 7) as follows:

BBI	HCVrep1b/Ava.1
BBII	HCVrep1b/Ava.5
BBIII	HCVrep1b/Huh.2
BBIV	HCVrep1b/Ava.7
BBV	HCVrep1b/Ava.2
BBVI	HCVrep1b/Clone A
BBVII	HCVrep1b/Clone B

5. Because the cell colonies were G418 resistant, we expected that the resistance was conferred by HCV replicons comprising adaptive mutations, harbored by those colonies. We tested this theory by sequencing the replicons, which were amplified from cDNA reverse transcribed from RNA isolated from each of the independent G418 resistant cell clones, before March 31, 2000. That data is presented at A3-A19. We then engineered each mutation back into the HCVrep1bBartMan/Avall backbone, as described in the Example in the Patent Application. We then transcribed RNA from each reconstructed replicon and electroporated it into naïve Huh7 cells, and compared the number of G418 resistant colonies compared to that obtained for the HCVrep1bBartMan/Avall replicon containing wild type NS5A (see A1, for example, where it was determined that the mutation identified in clone BBI was capable of increasing the frequency of G418 resistant colonies). Based on that result, we reasonably expected that the other mutations identified would similarly confer increased frequency of G418 resistant colonies, due to increased transfection efficiency of the mutant HCV.

7. I understand that willful false statements and the like are punishable by fine or imprisonment, or both (18 U.S.C. 1001) and may jeopardize the validity of the application or any patent issuing thereon.

Dr. Keril J. Blight

August __, 2004

Electroporation Huh7B p60

1 μ g vs 10 μ g — WT/Ava.I
 + 9 μ g — pol⁺/Ava.I
 Huh7B
 cellular — Original (HCVrepBarMan/Ava.I)
 RNA(4-1-00) — 10 μ g cellular RNA only

- * 10 T175's split 1:2 ~ 24hr prior to electroporation
- * Procedure as always (Plated 1/3 & 2/3 on p150's. Also removed 0.5ml of G18 for p60)
- * 4.2×10^7 cells total — Resuspend 3ml D-PBS $\Rightarrow 5.5 \times 10^6$ cells/eq
- * 24hr post-electroporation add G418 at 1mg/ml
- * Trypsinise p60 dishes & seed on 8 well chamber slides for IF (no G418)
- * Acetone fix ~ 8hr post-seeding
 - * Plated remainder of cells (mg/ml) on a p120. Added G418 ~ 16hr after seeding
- \Rightarrow IF for NS3 (H7) \Rightarrow NEGATIVE!

\Rightarrow The deletion is adaptive

Electroporation Huh7 (CMR) p49

1 μ g vs 10 μ g — WT/Ava.I
 + 9 μ g — pol⁺/Ava.I
 Huh7B
 cellular — Original (HCVrepBarMan/Ava.I)
 RNA(4-1-00) — 10 μ g cellular RNA only

- * 8 T175's split 1:2 25hr prior to electroporation
- * Procedure as previously (Plated 1/3 & 2/3 on p150's. Also removed 0.5ml from 1ml for p60)
- * 9.2×10^7 cells total \Rightarrow Resuspend 6ml PBS $\Rightarrow 6 \times 10^6$ cells/EP
- * 24hr post-electroporation add G418 at 0.8mg/ml

\Rightarrow Deletion is adaptive in CMR Huh7's. In fact, more colonies are consistently observed for HCVrepBarMan & HCVrep/Ava.I in CMR Huh7 cells vs Bortnischlager's.

Big Dye SeqClone A
~40ng PCR product

KB 486	A.1	CMR #829	2 μ l
KB 487	A.1	# 862	
488	A.1	# 869	
489	A.2	# 884	
490		# 885	
491		# 1038	
492		# 1039	
493	A.15	# 1039	1 μ l
494		# 1038	
495		# 949	
496		# 950	
497		# 970	
498		# 982	
499		# 971	
500		# 1030	
501	A.12	# 1042	1 μ l
502		# 931	
503		# 1047	
504		# 936	
505		# 923	
506		# 8990	
507	A.10	# 1040	1.2 μ l
508		# 1047	
509		# 936	
510		# 931	
511		# 923	
512	A.13	# 1029	4 μ l

KB 574	CMR #983	PCR # 15
575	# 1019	A.15
576	# 1022	A.15
577	# 1023	A.15
578	# 819	A.20
579	# 1030	A.15

Reamplified the gel
purified product
14-2-00
except KB578-A.2

Tuesday,

Untitled

Construction parameters:


Match Size	12
Maximum Added Gap Length in Contig	70
Maximum Added Gap Length in Sequence	70
Minimum Match Percentage	65
Maximum Register Shift Difference	70
Last group Considered	2
Gap Penalty	0.00
Gap Length Penalty	0.70
Consensus Threshold	75

Clone A

CREATING NEW contig 1: from Hartman WT(1>8012)
ENTERING 01-KB486(12>691) in contig 1: percent match 99
ENTERING 02-KB487(33>650) in contig 1: percent match 94
ENTERING 03-KB488(36>444) in contig 1: percent match 99
ENTERING 04-KB489(18>657) in contig 1: percent match 98
ENTERING 05-KB490(26>629) in contig 1: percent match 95
ENTERING 06-KB491(1>516) in contig 1: percent match 98
ENTERING 07-KB492(8>848) in contig 1: percent match 95
ENTERING 08-KB493(1>520) in contig 1: percent match 97
ENTERING 09-KB494(2>740) in contig 1: percent match 97
ENTERING 10-KB495(1>742) in contig 1: percent match 96
ENTERING 11-KB496(21>786) in contig 1: percent match 98
ENTERING 12-KB497(20>737) in contig 1: percent match 96
ENTERING 13-KB498(4>843) in contig 1: percent match 96
ENTERING 14-KB499(1>661) in contig 1: percent match 93
ENTERING 15-KB500(11>779) in contig 1: percent match 97
ENTERING 16-KB501(1>112) in contig 1: percent match 92
ENTERING 17-KB502(17>791) in contig 1: percent match 98
ENTERING 18-KB503(1>752) in contig 1: percent match 99
ENTERING 19-KB504(2>737) in contig 1: percent match 97
ENTERING 20-KB505(21>720) in contig 1: percent match 98
ENTERING 21-KB506(8>757) in contig 1: percent match 98
ENTERING 22-KB507(1>633) in contig 1: percent match 98
ENTERING 23-KB508(1>728) in contig 1: percent match 98
ENTERING 24-KB509(7>745) in contig 1: percent match 98
ENTERING 25-KB510(18>703) in contig 1: percent match 99
ENTERING 26-KB511(1>711) in contig 1: percent match 94
Sequence 27-KB512 was not added, it is all poor data
Elapsed Time 0:0:16

- * Separate PCR products on 0.8% TBE
- * Isolate DNA from gel slices (as on 5-2-00)

- * Pellet purified DNA & wash 2x 80% EtOH
- * Resuspend 10µl TE
- * Quantitate 0.5µl



Big Dye Seq ~ 100ng PCR product / 0.8µl oligo @ 4pmol/µl

KB513	CMR #862	PCR #1	KB528	CMR #1029	PCR #3	KB571	CMR
514	#867		529	#1022	PCR #3	572	CMR
515	#829		530	#1030	PCR #3	573	CMR
516	#884	PCR #2	531	#1040	PCR #4		
517	#885		532	#1047			
518	#1038		533	#936			
519	#1039		534	#931			
520	#1038	PCR #3	535	#923			
521	#1039		536	#1047	PCR #5		
522	#949		537	#936			
523	#950		538	#931			
524	#910		539	#923			
525	#971		540	#1042			
526	#482		541	#899c			

Wednesday,

Untitled

Construction parameters:

Match Size	12
Minimum Added Gap Length in Contig	70
Minimum Added Gap Length in Sequence	70
Minimum Match Percentage	65
Maximum Register Shift Difference	70
Lastgroup Considered	2
Gap Penalty	0.00
Gap Length Penalty	0.70
Consensus Threshold	75

Avo.1

CREATING NEW contig 1: from BartMan WT copy(1>8012)

ENTERING 01.KB513(16>572) in contig 1: percent match 96

ENTERING 02.KB514(33>443) in contig 1: percent match 99

ENTERING 03.KB515(1>600) in contig 1: percent match 99

ENTERING 04.KB516(1>806) in contig 1: percent match 95

ENTERING 05.KB517(5>698) in contig 1: percent match 94

ENTERING 06.KB518(1>518) in contig 1: percent match 98

ENTERING 07.KB519(14>821) in contig 1: percent match 94

ENTERING 08.KB520(6>711) in contig 1: percent match 96

ENTERING 09.KB521(11>514) in contig 1: percent match 98

ENTERING 10.KB522(14>666) in contig 1: percent match 97

ENTERING 11.KB523(11>756) in contig 1: percent match 97

ENTERING 12.KB524(22>729) in contig 1: percent match 95

ENTERING 13.KB525(24>734) in contig 1: percent match 92

ENTERING 14.KB526(1>707) in contig 1: percent match 97

ENTERING 15.KB527(19>853) in contig 1: percent match 91

ENTERING 16.KB528(23>535) in contig 1: percent match 85

ENTERING 17.KB529(43>658) in contig 1: percent match 97

ENTERING 18.KB530(42>775) in contig 1: percent match 89

ENTERING 19.KB531(7>679) in contig 1: percent match 98

ENTERING 20.KB532(19>723) in contig 1: percent match 98

ENTERING 54.KB533(1>859) in contig 1: percent match 96

ENTERING 55.KB534(1>750) in contig 1: percent match 97

ENTERING 56.KB535(17>755) in contig 1: percent match 98

ENTERING 57.KB536(8>829) in contig 1: percent match 97

ENTERING 58.KB537(1>763) in contig 1: percent match 96

ENTERING 59.KB538(12>697) in contig 1: percent match 99

ENTERING 60.KB539(18>739) in contig 1: percent match 99

ENTERING 61.KB540(22>114) in contig 1: percent match 89

ENTERING 62.KB541(25>692) in contig 1: percent match 98

Elapsed Time 0:0:18

8.2 8.3 A.15 100ng DNA

B.2 80ng/μl

B.3 80ng/μl

A.15 120ng/μl

Big Dye Seq - 40ng PCR product

KB542	CMR #862	B.1	KB560	CMR #1040	B.4
543	#867	B.1	561	#1047	
544	#879	B.1	562	#936	
545	#884	B.2	563	#931	
546	#885		564	#923	
547	#1038		565	#1047	B.5
548	#1039		566	#936	
549	#1038	B.3	567	#931	
550	#1039		568	#923	
551	#949		569	#1042	
552	#950		570	#999c	
553	#970				
554	#971				
555	#982				
556	#983				
557	#1029				
558	#1022				
559	#1030				

Thursday,
Untitled

Construction parameters:	12
Match Size	70
Maximum Added Gap Length in Contig	70
Maximum Added Gap Length in Sequence	65
Minimum Match Percentage	70
Maximum Register Shift Difference	2
Last group Considered	0.00
Gap Penalty	0.70
Gap Length Penalty	75
Consensus Threshold	

Clone B

CREATING NEW contig 1: from BartMan WT copy(1>8012)

ENTERING 31-KB542(41>556) in contig 1: percent match 99

ENTERING 32-KB543(34>444) in contig 1: percent match 99

ENTERING 33-KB544(3>597) in contig 1: percent match 99

ENTERING 34-KB545(1>690) in contig 1: percent match 95

ENTERING 35-KB546(1>598) in contig 1: percent match 96

ENTERING 36-KB547(1>515) in contig 1: percent match 97

ENTERING 37-KB548(1>574) in contig 1: percent match 95

ENTERING 38-KB549(12>678) in contig 1: percent match 98

ENTERING 39-KB550(29>512) in contig 1: percent match 98

ENTERING 40-KB551(15>620) in contig 1: percent match 97

ENTERING 41-KB552(7>676) in contig 1: percent match 98

ENTERING 42-KB553(12>594) in contig 1: percent match 98

ENTERING 43-KB554(2>597) in contig 1: percent match 93

ENTERING 44-KB555(1>696) in contig 1: percent match 98

ENTERING 45-KB556(5>707) in contig 1: percent match 96

ENTERING 46-KB557(6>567) in contig 1: percent match 95

ENTERING 47-KB558(37>700) in contig 1: percent match 96

ENTERING 48-KB559(18>596) in contig 1: percent match 96

ENTERING 49-KB560(2>646) in contig 1: percent match 97

ENTERING 50-KB561(7>662) in contig 1: percent match 97

ENTERING 51-KB562(1>677) in contig 1: percent match 98

ENTERING 52-KB563(1>693) in contig 1: percent match 97

ENTERING 53-KB564(16>725) in contig 1: percent match 97

ENTERING 54-KB565(1>695) in contig 1: percent match 98

ENTERING 55-KB566(8>687) in contig 1: percent match 98

ENTERING 56-KB567(1>600) in contig 1: percent match 99

ENTERING 57-KB568(15>732) in contig 1: percent match 97

ENTERING 58-KB569(14>127) in contig 1: percent match 85

ENTERING 59-KB570(7>686) in contig 1: percent match 98

Elapsed Time 0:0:17



Page 10

PAGE 21/37 * RCVD AT 8/30/2004 6:12:40 PM [Eastern Daylight Time] * SVR:USPTO-EFXRF-1/5 * DNIS:8729308 * CSID: * DURATION (mm-ss):11-12

Big Dye Seq ~40ng PCR product Ava.5

Ava.5 PCR #1 30ng/μl
 PCR #2 170ng/μl
 PCR #3 150ng/μl
 PCR #4 60ng/μl
 PCR #5 30ng/μl

iscalculated
 Added
 5ng PCR
 adbot, instead
 of 40ng

KB580	CMR #862	PCR #1
581	#867	↓
582	#829	↓
583	CMR #884	PCR #2
584	#885	↓
585	#1038	↓
586	#1039	↓
587	#1014	↓
588	CMR #1038	PCR #3
589	#1039	↓
590	#949	↓
591	#950	↓
592	#970	↓
593	#971	↓
594	#982	↓
595	#983	↓
596	#1029	↓
597	#1022	↓
598	#1030	↓
599	#1023	↓

KB600	CMR #1040	PCR #4
601	#1047	↓
602	#936	↓
603	#931	↓
604	#923	↓
605	#1047	PCR #5
606	#936	↓
607	#931	↓
608	#913	↓
609	#1042	↓
610	#849C	↓

Tuesday,

Untitled

Construction parameters:	12
Match Size	70
Minimum Added Gap Length in Contig	70
Minimum Added Gap Length in Sequence	65
Minimum Match Percentage	70
Maximum Register Shift Difference	2
Lastgroup Considered	0.00
Gap Penalty	0.70
Gap Length Penalty	75
Consensus Threshold	

Ava.5

CREATING NEW contig 1: from PartMan Wf copy 1(1>8012)

ENTERING 01.KB580(33>533) in contig 1: percent match 97

ENTERING 02.KB581(1>439) in contig 1: percent match 94

ENTERING 03.KB582(16>598) in contig 1: percent match 90

Sequence 04.KB583 was not added, it is all poor data

Sequence 05.KB584 was not added, it is all poor data

Sequence 06.KB585 was not added, it is all poor data

NOT ENTERING in contig 1: 07.KB586(11>408) due to percent match (43) below threshold 65

NOT ENTERING in contig 1: 07.KB586(11>488) due to percent match (47) below threshold 65

CREATING NEW contig 2: from 07.KB586(11>488)

Sequence 08.KB587 was not added, it is all poor data

ENTERING 09.KB588(1>665) in contig 1: percent match 97

ENTERING 10.KB589(8>515) in contig 1: percent match 98

ENTERING 11.KB590(1>739) in contig 1: percent match 96

ENTERING 12.KB591(1>734) in contig 1: percent match 96

ENTERING 13.KB592(8>738) in contig 1: percent match 95

ENTERING 14.KB593(1>576) in contig 1: percent match 94

ENTERING 15.KB594(1>708) in contig 1: percent match 95

ENTERING 16.KB595(9>698) in contig 1: percent match 94

ENTERING 17.KB596(3>558) in contig 1: percent match 94

ENTERING 18.KB597(34>531) in contig 1: percent match 98

ENTERING 19.KB598(20>545) in contig 1: percent match 98

ENTERING 20.KB599(5>325) in contig 1: percent match 94

ENTERING 21.KB600(1>596) in contig 1: percent match 97

NOT ENTERING in contig 2: 22.KB601(1>701) due to percent match (51) below threshold 65

ENTERING 22.KB601(1>701) in contig 1: percent match 96

NOT ENTERING in contig 2: 23.KB602(1>540) due to percent match (50) below threshold 65

ENTERING 23.KB602(1>540) in contig 1: percent match 98

NOT ENTERING in contig 2: 24.KB603(6>648) due to percent match (45) below threshold 65

ENTERING 24.KB603(6>648) in contig 1: percent match 98

ENTERING 25.KB604(13>733) in contig 1: percent match 95

NOT ENTERING in contig 2: 26.KB605(1>735) due to percent match (51) below threshold 65

ENTERING 26.KB605(1>735) in contig 1: percent match 98

NOT ENTERING in contig 2: 27.KB606(1>603) due to percent match (50) below threshold 65

ENTERING 27.KB606(1>603) in contig 1: percent match 99

NOT ENTERING in contig 2: 28.KB607(1>703) due to percent match (45) below threshold 65

ENTERING 28.KB607(1>703) in contig 1: percent match 95

ENTERING 29.KB608(17>727) in contig 1: percent match 96

ENTERING 30.KB609(14>117) in contig 1: percent match 90

ENTERING 31.KB610(43>698) in contig 1: percent match 97

Elapsed Time 0:0:17



Page 5

[illegible][illegible][illegible][illegible][illegible]

Big Dye Seq ~ 40ng PCR product Ava.2

Ava.2 PCR #1 20ng/ μ l
 PCR #2 120ng/ μ l
 PCR #3 80ng/ μ l
 PCR #4 40ng/ μ l
 PCR #5 70ng/ μ l

	CMR #862	PCR #1	KB631	CMR #1040	PCR #4
KB611	#867	↓	KB632	#1041	↓
612	#829	↓	KB633	#936	↓
613	#884	PCR #2	KB634	#931	↓
* 614	#885	↓	KB635	#923	↓
* 615	#1038	↓	KB636	#1047	PCR #5
* 616	#1039	↓	637	#936	PCR #5
617	#1014	↓	638	#931	↓
* 618	#1038	PCR #3	639	#923	↓
619	#1039	↓	640	#1042	↓
620	#949	↓	641	#899c	↓
621	#950	↓			
622	#770	↓			
623	#971	↓			
624	#982	↓			
625	#983	↓			
626	#1029	↓			
627	#1022	↓			
628	#1030	↓			
629	#1023	↓			
630		↓			

* Reactions 614-616 & 618 did not work → See Emily's comments



Wednesday,

Project: Unified Contlg 1

Page 9

Barthman WT copy (1>8012)	5180	5190	5200	5210	5220	5230	5240	5250	5260
18-KB628 (11>568)									
19-KB629 (11>702)									
17-KB527 (8>568)									
20-KB630 (23>335)									
Barthman WT copy (1>8012)	5270	5280	5290	5300	5310	5320	5330	5340	5350
19-KB629 (11>702)									
17-KB527 (8>568)									
20-KB630 (23>335)									
Barthman WT copy (1>8012)	5360	5370	5380	5390	5400	5410	5420	5430	5440
19-KB629 (11>702)									
17-KB527 (8>568)									
20-KB630 (23>335)									
Barthman WT copy (1>8012)	5460	5470	5480	5490	5500	5510	5520	5530	5540
19-KB629 (11>702)									
17-KB527 (8>568)									
20-KB630 (23>335)									
Barthman WT copy (1>8012)	5550	5560	5570	5580	5590	5600	5610	5620	5630
19-KB629 (11>702)									
17-KB527 (8>568)									
20-KB630 (23>335)									
21-KB631 (2>679)									
Barthman WT copy (1>8012)	5650	5660	5670	5680	5690	5700	5710	5720	5730
21-KB631 (2>679)									

Wednesday.

Untitled

Construction parameters:

Vector Size	12
Minimum Added Gap Length in Contig	70
Minimum Added Gap Length in Sequence	70
Minimum Match Percentage	65
Maximum Register Shift Difference	70
Lastgroup Considered	2
Gap Penalty	0.00
Gap Length Penalty	0.70
Consensus Threshold	75

Ava. 2

CREATING NEW contig 1: from BartMan WT copy(1>8012)

ENTERING 01-KB611(36>680) in contig 1: percent match 95

ENTERING 02-KB612(35>446) in contig 1: percent match 99

ENTERING 03-KB613(11>666) in contig 1: percent match 99

Sequence 04-KB614 was not added, it is all poor data —

Sequence 05-KB615 was not added, it is all poor data —

Sequence 06-KB616 was not added, it is all poor data —

NOT ENTERING in contig 1: 07-KB617(1>597) due to percent match (47) below threshold 65

NOT ENTERING in contig 1: 07-KB617(1>597) due to percent match (51) below threshold 65

CREATING NEW contig 2: from 07-KB617(1>597)

Sequence 08-KB618 was not added, it is all poor data

NOT ENTERING in contig 2: 09-KB619(16>662) due to percent match (45) below threshold 65

ENTERING 09-KB619(16>662) in contig 1: percent match 92

ENTERING 10-KB620(9>517) in contig 1: percent match 96

ENTERING 11-KB621(65>627) in contig 1: percent match 96

Sequence 12-KB622 was not added, it is all poor data —

ENTERING 13-KB623(38>575) in contig 1: percent match 97

ENTERING 14-KB624(6>414) in contig 1: percent match 92

ENTERING 15-KB625(1>700) in contig 1: percent match 98

ENTERING 16-KB626(20>542) in contig 1: percent match 95

ENTERING 17-KB627(8>568) in contig 1: percent match 94

ENTERING 18-KB628(38>696) in contig 1: percent match 96

ENTERING 19-KB629(11>782) in contig 1: percent match 96

ENTERING 20-KB630(23>335) in contig 1: percent match 95

ENTERING 21-KB631(2>679) in contig 1: percent match 98

NOT ENTERING in contig 2: 22-KB632(1>707) due to percent match (50) below threshold 65

ENTERING 22-KB632(1>707) in contig 1: percent match 97

NOT ENTERING in contig 2: 23-KB633(1>721) due to percent match (51) below threshold 65

ENTERING 23-KB633(1>721) in contig 1: percent match 97

NOT ENTERING in contig 2: 24-KB634(7>642) due to percent match (47) below threshold 65

ENTERING 24-KB634(7>642) in contig 1: percent match 98

ENTERING 25-KB635(13>726) in contig 1: percent match 97

NOT ENTERING in contig 2: 26-KB636(1>680) due to percent match (50) below threshold 65

ENTERING 26-KB636(1>680) in contig 1: percent match 98

NOT ENTERING in contig 2: 27-KB637(1>688) due to percent match (50) below threshold 65

ENTERING 27-KB637(1>688) in contig 1: percent match 97

NOT ENTERING in contig 2: 28-KB638(12>658) due to percent match (45) below threshold 65

ENTERING 28-KB638(12>658) in contig 1: percent match 98

ENTERING 29-KB639(15>728) in contig 1: percent match 97

ENTERING 30-KB640(8>135) in contig 1: percent match 85

ENTERING 31-KB641(1>666) in contig 1: percent match 86

Elapsed Time 0:0:16

DNA Strider™ 1.317 ### Wednesday,

04:30 PM

Clone A n.t. 5336

II Clone B n.t. 5336

I Ant 5345-5485

II Ava. 2 n.t. 5320

II Ava. 5 n.t. 3550 & n.t. 4573
n.t. 5290

HCVrep1b BartMan/AvaII [1801 to 7758] -> Translate - 1-frame

DNA sequence 11313 bp gccaggaacccga .7. cgactccatata circular

artemischlager replicon I377/NS3-3'UTR (Genbank AJ242652).

Constructed in the pMT backbone.

Marked by AvaII in the variable region of the 3'UTR.

NS3

1801/1
atg ggc cct att acg gcc tac tcc caa cag acg cga ggc cta ctt ggc tgc atc atc act
M A V I T A Y S Q Q T R Q L L G Q I I T

1861/21
agg ctc aca ggc cgg gac agg aac cag gtc gag ggg gag gtc caa gty gtc tcc acc gca
S L T Q R D R H Q V R G E V Q V V S T A

1921/41
aca caa tct ttc ctg ggc acc tgc gtc aat ggc gtc tgt tgg act gtc tat cat ggt gcc
T Q S F L A T Q V N G V C W T V Y H Q A

1981/61
ggc tca aag acc ctt gcc ggc cca aag ggc cca atc acc caa atg tac aac aat gtc gac
G S K T L A G E K G P I T Q M T T N V D

2041/81
cag gaa ctc gtc ggc tgc ggc cca ggc ggc ggc ggc ggc ggc ggc ggc ggc ggc ggc ggc
Q D L V G W Q A P P G A R S L R F C T C

2101/101
ggc agt tgc aac ctt tac tgc ggc acc ggc ggc ggc ggc ggc ggc ggc ggc ggc ggc ggc
G S S D D T L V T N H A D V I P V R R R

2161/121
ggc gaa acc ggc ggc ggc ggc ggc ggc ggc ggc ggc ggc ggc ggc ggc ggc ggc ggc
G D S R Q S L L S P R P V A Y L E Q S S

2221/141
ggc ggt cca ctg ctc tgc ccc tgc ggc cca acc ggc ggc ggc ggc ggc ggc ggc ggc ggc
G G D L L C P S Q M A V G I F R A A V C

2281/161
ggc ggc ggc ggc ggc ggc ggc ggc ggc ggc ggc ggc ggc ggc ggc ggc ggc ggc ggc
R G V A K A V D F V P V H S M E T G M

2341/181
ggc tcc cgc gtc ttc acc gaa aac cgc tcc ctt cgc ggc gta cgc acc ttc cag gtc
K S F V P T D M B S P P A V P Q T F Q V

2401/201
ggc cat cta caa gcc cct act ggt agc ggc aag agc act aag gtc ccc ggt ggt tat gca
A H L H A P P Q S G K S T K V P A A T A

2461/221
ggc cca ggc tat aag gtc ctt gtc ctg aac cgc tcc gtc gcc gcc acc cta ggt ttc ggc
A Q G Y K V L V L N P S V A A T L Q F G

2521/241
ggc tat atg tct aag gca cat ggt acc gac cct aac atc aga acc ggc gta agc acc atc
A Y M S K A R Q I D P M I R T G V R T I

2581/261
acc acc ggt gcc ccc atc acc tac tcc acc tat ggc aag ttt ctt gcc gac ggt ggt tgc
S T Q A P I T T S T Y G K P L A D G G C

2641/281
tct ggc ggc gcc tat gac atc ata ata tgc gat gag tgc cac tca acc gac tgc acc act
S G Q A T D I I I C D E C H S T D S T T

2701/301
atc ctg ggc atc ggc aca gtc ctg gac caa ggc gag acc ggt gga ggc cga ctc gtc gtc
I L G I G T V L D Q A E T A G A R L V V

2761/321
ctc gcc acc ggt acc cct cgc ggc tgc gtc acc gtc cca cat ccc aag atc gag gag gtc
L A T A T P P Q S V T V P H P H I E E V

2821/341
ggt ctg tcc agc act gga gaa atc ccc ttt tat ggc aag gcc atc ccc atc gaa acc atc
A L S S T G H I P P Y G K A I P I E T I

2881/361
ggc ggc ggc ggc ggc ggc ggc ggc ggc ggc ggc ggc ggc ggc ggc ggc ggc ggc ggc
G G G R H L I P C H S K K K C D E L A A

2941/381
aag ctg tcc ggc ctc gga ctc aat ggt gta gca tat tac cgc ggc ctt gat gta tcc gtc
K L S G L G L H A V A Y Y R G L D V S V

3001/401
atc cca act aac gga gac gcc att ttc gta gca aag gac ggc cta atg acc ggc ttt acc
T T C C A A G G A G C A T T T C G A G C A G G C T A T G A C G G C T T A C

Ava. II/Clone C/Clone D
n.t. 5336

Ava. 13 n.t. 5320

Ava. 7 n.t. 5313

Huh. 2 n.t. 5314

Clone A
cga → ggc 87
Gln ArgE → G 177
GAA GGAT → I 255
ACC ATC

CloneAva.5
ag⁺ → gg⁺
Ser Gly

Huh. 5.

$O = 958 \text{ sites}$
1172/1177/1179

Clone Hh. 4
Clone Ava. 11
Clone Ccd
Clone A. 4 B
agg → atc
Ser Ile

Clone Ava
agg → gg
Arg Gl

gcc → tcc
Ala Ser
Clone Ava. 2
Ava. 13

447 aa Clone Ave. I
n. + 5345-5425

ISDR (4000)

① Amino acid changes
of 90 con

$TCR \rightarrow CEC$
 Ser Pro

1447.52

1142

19-3-00

G418-colonies picked for BartMan/AvaII in CMR Huh7 cells

Huh.4 10 μ g. $\frac{1}{3}$ 8-2-00
 Huh.5 1 μ g $\frac{2}{3}$ "
 Huh.6 1 μ g $\frac{1}{3}$ "
 Huh.7 1 μ g $\frac{1}{3}$ "
 Huh.8 10 μ g $\frac{1}{3}$ "
 Huh.9 10 μ g $\frac{1}{3}$ "

⇒ G418 added to 750 μ g/ml on 20-3-00Huh.4

- Transfer to 24 well plate 21-3-00 p1
- Transfer to 6 well plate 22-3-00 p2
- Transfer to T25 flask 25-3-00 p3

ALSO, Cell count #1 3.2×10^5 cells/ml } $\sim 3 \times 10^5$ cells/ml
 Cell count #2 2.95×10^5 "

⇒ Trizol extract RNA from 80 μ l cells (2.4×10^4 cells)

- Transfer to T75 flask 30-3-00 p4

10-4-00 - Froze 12 vials @ 4×10^6 cells/vial Huh.4 p7 (split 1:2 8-4-00)

Tank 2 Rack 5 Box 1 (1 vial); Box 6 (1 vial); Box 7 (1 vial);

Box 8 (6 vials); Box 9 (1 vial) & Tank 2 Rack 7 Box 3 (2 vials)

Huh.5

- Transfer to 24 well plate 24-3-00 p1
- Transfer to 12 well plate 28-3-00 p2
- Transfer to 6 well plate 30-3-00 p3
- Transfer to T25 flask 2-4-00 p4

ALSO, Cell count #1 4.4×10^5 cells/ml } 4.3×10^5
 Cell count #2 4.15×10^5 "

⇒ Trizol extract RNA from 80 μ l ($\sim 3.4 \times 10^4$ cells)

- Transfer to T75 flask 2-4-00 p5

20-4-00 - Froze 10 vials @ 4×10^6 cells/vial p8 (3 T75s split 1:1.5 on 11-4-00)

Huh. 7

- Transfer to 24 well plate 22-3-00 p1
- Transfer to 12 well plate 25-3-00 p2
- Transfer to 6 well plate 26-3-00 p3
- Transfer to T25 flask 29-3-00 p4

ALSO, Cell count #1 2.75×10^5 cells/ml } 3×10^5 cells/ml
 Cell count #2 3.5×10^5 }

⇒ Trizol extract RNA from 80µl (2.4×10^4 cells)

- Transfer to T75 flask 2-4-00 p5

12-4-00 - Froze 12 vials @ 4×10^6 cells/vial Huh. 7 p8 (split 1:2 on 10-4-00)

Tank 2 Rack 7 Box 3 (2 vials)

" " Box 4 (7 vials)

" " Box 5 (3 vials)

Huh. 8

- Transfer to 24 well plate 21-3-00 p1
- Transfer to 12 well plate 24-3-00 p2
- Transfer to 6 well plate 26-3-00 p3
- Transfer to T25 flask 29-3-00 p4

ALSO, Cell count #1 1.4×10^5 cells/ml } 1.75×10^5 cells/ml
 Cell count #2 2×10^5 }

⇒ Trizol extract RNA from 85µl (1.5×10^4 cells)

- Transfer to T75 flask 2-4-00 p5

14-4-00 Froze 10 vials @ 4×10^6 cells/vial Huh. 8 p8 (split 1:2 12-4-00)

Tank 2 Rack 3 Box 2 (1 vial)

Box 4 (1 vial)

Box 5 (1 vial)

Box 6 (1 vial)

Tank 1 Rack 3 Box 4 (1 vial)

Box 8 (1 vial)

Remainder stored in -80°C freezer box

Hub.9

- Transfer to 24 well plate 28-3-00 p1
 - Transfer to 12 well plate 6-4-00 p2 (media change 12-4-00)

- Transfer to 6 well plate 14-1-00 p3

- Transfer to T25 flask 19-4-00 p4

ALSO, Cell count #1 6.9×10^5 cells/ml } $\sim 8.5 \times 10^5$ cells/ml
 Cell count #2 10^6 cells/ml }

→ Trizol extract RNA from 70 μ l ($\sim 6 \times 10^4$ cells)

- Transfer to T75 flask 23-4-00 p5

10-5-00 - Freeze 12 vials at 4×10^6 cells/vial p8

Stored -80°C freezer box

Clone D

- Transfer to 24 well plate p1

- Transfer to 12 well plate p2

- Transfer to 6 well plate p3

- Transfer to T25 flask 16-3-00 p4

ALSO, Cell count #1 3.4×10^5 cells/ml } 3×10^5 cells/ml
 Cell count #2 2.5×10^5 }

→ Trizol extract 80 μ l (2.7×10^4 cells)

- Transfer to T75 flask 20-3-00 p5

31-3-00 - Freeze 7 vials @ 3.5×10^6 cells/vial p8 (4 T75's split 1:2 on 29-3-00)

Tank 2 Rack 7 Box 7

20-4-00Electroporation Huh7(cmr) p69~ 6×10^6 cells/electroporation+ 1.5 μ g replicon RNA + 5 μ g Huh7B cellular RNA (4-1-00)

- (1) HCVrep/Ava.2 BB.V
- (2) HCVrep/Ava.5 BB.II
- (3) HCVrep/Clone B BB.VII
- (4) HCVrep/Ava.1 BB.I
- (5) HCVrep13/ Δ ISDR
- (6) HCVrep BartMao/Ava.II (original)
- (7) HCVrep(pol-)/Ava.1
- (8) Polio rep-GFP

Place electroporated cells into media. Total volume = 9.5ml

A. Plate 0.5ml / p100

B. Plate 3ml & 6ml per p150

21-4-00

* At G418 at 0.8mg/ml ~ 27hr post-electroporation

* Stain p100's on 6-5-00


* Stain p100 with 3ml plated on 14-5-00

N.B. When the cells are plated too dense they begin to detach due to overconfluence, particularly Clone B & Ava. 2 which appear to replicate in >90% of Huh7 cells


* Clone B > Ava. 2 > Ava. 5 & Ava. 1 # of G418-resistant colonies

* No colonies observed for H77 HCVrep13/ Δ ISDR & HCVrep(pol-)/Ava. 1


G418-resistant colonies (Experiment 20-4-00)




II




HCVrep13/AISDR




V



pol-



VII




HCVrep1bBartMan/Avall



I

Check 1µg (~2µg) on non-denaturing gel



21-5-00

Electroporation Huh7 p69

5 T175's split 1:2 ~ 24hr prior

1.5µg HCVrep RNA + 4µg Huh7B cellular RNA (4-1-00)

~ 6×10^6 cells/electroporation

1. HCVrep1b/Ava.1 BB I
2. HCVrep1b/Ava.2 BB V
3. HCVrep1b/Ava.5 BB II
4. HCVrep1b/Ava.7 BB IV
5. HCVrep1b/Clone B BB VII
6. HCVrep1b/Huh.2 BB III
7. HCVrep1b/Ser→Ile (1179)
8. HCVrep1bBartMan/AvaII (original)
9. HCVrep1bBartMan(pol⁻)/AvaII

Volume total ~ 9.5ml → Plate 0.4ml, 0.3ml, 0.2ml, 0.1ml per p100

Plate 1.5ml per p150

(N.B. Also plated 6ml/p150 for HCVrep1b/Ser→Ile)

22-5-00Electroporation HeLa cells p15

* 5 T175's split 30hr prior 1:2

* 1.3×10^8 cells total resuspended in 2.5ml ice-cold D-PBS* 99 μ sec, 0.9 kV, 5 pulses, 0.4ml cells ($\sim 6 \times 10^6$ cells)2 μ g HCVrep RNA

1 HCVrep1b/Clone B BB VII

2 HCVrep1b/Ava.1 BB I

3 HCVrep1b/Ava.2 BB V

4 HCVrep1b/Ava.5 BB II

5 HCVrep1b/Ava.7 BB IV

6 HCVrep1b/Huh.2 BB III

7 HCVrep1b/BartMan/AvaII

8 HCVrep1b/BartMan(pol-)/AvaII

9 HCVrep13/S \rightarrow I10 HCVrep13/ Δ ISDR

11 No RNA

* $V_T \sim 9.5$ ml media + EPed cells

* Plate 3ml & 6ml per p150

" 0.5ml per p100 (For no RNA EP, plated total cells on p150)

24-5-00

At 48hr post-EP add 0.9 mg/ml G418

 \Rightarrow No colonies observed for any of the RNAs electroporated

**This Page is Inserted by IFW Indexing and Scanning
Operations and is not part of the Official Record**

BEST AVAILABLE IMAGES

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images include but are not limited to the items checked:

- ☐ **BLACK BORDERS**
- ☐ **IMAGE CUT OFF AT TOP, BOTTOM OR SIDES**
- ☒ **FADED TEXT OR DRAWING**
- ☒ **BLURRED OR ILLEGIBLE TEXT OR DRAWING**
- ☐ **SKEWED/SLANTED IMAGES**
- ☐ **COLOR OR BLACK AND WHITE PHOTOGRAPHS**
- ☐ **GRAY SCALE DOCUMENTS**
- ☐ **LINES OR MARKS ON ORIGINAL DOCUMENT**
- ☐ **REFERENCE(S) OR EXHIBIT(S) SUBMITTED ARE POOR QUALITY**
- ☐ **OTHER:** _____

IMAGES ARE BEST AVAILABLE COPY.

As rescanning these documents will not correct the image problems checked, please do not report these problems to the IFW Image Problem Mailbox.